

Pulson 100

The Pulson 100, for ultrasound therapy.

The Pulson : a 'must' for your practice. Portable and supremely easy to use, with a large screen. Pushbuttons for the most important functions.



Article number: 332.343

Characteristics

- Ultrasound :
 - 1 Frequency head (1 MHz), 4 cm²
 - Detachable head
 - Pulsed and continuous mode (10 - 20 - 30 - 40 - 50 - 100%)
 - Acoustic and visual contact control
- Functionalities :
 - Indication: 11
 - Free mode

Technical specifications

Mains voltages	: 100- 240-VAC, 50/60 Hz +/- 10 %
Maximum power -in operation	: 35 VA
Dimensions (b x h x d)	: 266 x 100 x 275 mm
Weight including accessories	: ca 4,150 kg
Safety class	: Class I
Isolation type	: Type BF
Classification according to MDD	: Class IIa
Safety examination	: CE from conformity with Directives MDD 93/42/EEC

Standard accessories

- 100 689 Mains lead
- 332 310 US-head, 1 MHz, ERA 4 cm², incl. holder
- 100 016 Contact gel, 500 ml
- 117 277 Sets of photos US-head placement

Manuals

- 323 011 Safety Instructions
- 331 958 User Manual English
- 311 872 CD-ROM User Manuals multi language

Optional accessories

- 332 321 US-head, 1 MHz, ERA 1 cm²
- 302 955 Carrying case 200 series
- 100 019 Contact-gel, 5L
- 100 020 Pump, 5L
- 117 210 Mobil 700, trolley, 2 shelves, 1 drawer
- 117 220 Extra shelf + drawer for Mobil 700
- 117 215 Bottle support for Mobil 700