

Pulson 200

The Pulson 200, for ultrasound therapy.

The Pulson : a 'must' for your practice. Portable and supremely easy to use, with a large backlight screen. Pushbuttons for the most important functions.



Article number: 300 399

Characteristics

- Ultrasound :
 - Multifrequency head (1 and 3 MHz), 4 cm²
 - Continuous and pulsed mode (10 - 20 - 30 - 40 - 50 - 100%)
 - Acoustic contact control

- Functionalities :
 - Objectives: 14
 - Indication: 23
 - Diagnostic: 1
 - Free memory: 20

Technical specifications

Languages	: English, German, French, Dutch, Spanish, Italian, Portuguese and Polish
Mains voltages	: 100- 240-VAC, 50/60 Hz 50/60 HZ +/- 10 %
Maximum power -in operation	: 85 VA
Dimensions (b x h x d)	: 266 x 100 x 275 mm
Weight	: ca 3,650 kg
Safety class	: Class I
Isolationtype	: Type BF
Classification according to MDD	: Class IIa
Safety examination	: CE from conformity with Directives MDD 93/42/EEC

Standard accessories

- 100 689 Mains lead
- 117 122 US-head, multifrequent, ERA 4 cm², incl. holder
- 100 016 Contact gel, 500 ml
- 117 277 Sets of photos Pulson 200 US-head placement

Manuals

- 323 011 Safety Instructions
- 117 272 Treatment Guide Dutch
- 117 273 Treatment Guide French
- 117 275 Treatment Guide English
- 117 274 Treatment Guide German
- 327 393 Treatment Guide Spanish

Optional accessories

- 117 123 US-head, multi-freq., ERA 1 cm²
- 302 955 Carrying case 200 series
- 100 016 Contact gel, 500 ml
- 100 019 Contact gel, 5L
- 100 020 Pump, 5L
- 117 210 Mobil 700, trolley, 2 shelves, 1 drawer
- 117 220 Extra shelf + drawer for Mobil 700
- 117 215 Bottle support for Mobil 700