

Shock waves: a new type of physiotherapy treatment.

NEW



Storz Medical Masterpuls MP100

For the past 30 years, **GymnaUniphy** has been one of the largest manufacturers and suppliers of physical therapy and re-education devices. The company is the **exclusive worldwide distributor** of the Masterpuls MP 100, manufactured by Storz Medical, Switzerland's shock wave specialist. With the introduction of the **Storz Masterpuls MP100**, GymnaUniphy has now made shock wave therapy available to all physiotherapists.

Advantages of this Device

Shock wave therapy is a recent development that is also being used successfully by orthopaedic specialists. As compared to other treatment methods, the advantages for your patients are:

- Solution for chronic pathologies
- Avoid surgery
- Painless
- No medication needed
- No risk of allergy
- Fast, ambulant treatment
- Almost no side effects (bruise)
- Quickly reduces pain
- Only 3-5 treatments needed
- With the MP100, no medical imaging or anaesthesia is needed

Virtually painless treatment.



How does the Masterpuls MP100 work?

Shock waves are produced by pneumatic pulses that create ballistic movements. The shock waves are transmitted to the surrounding tissues. Shock waves are applied via a light hand piece that is easy to use, which is placed directly in contact with the skin at the site of the lesion. The effect is mainly mechanical.

For patients, the Masterpuls MP100 provides **fast, reliable** and **effective** treatment.

The Masterpuls MP100 can also be an **alternative to surgery**, particularly for patients suffering from chronic pain. The following pathologies can be successfully treated:

- Bursitis
- Tendinitis (golf and tennis elbow, patella tendon syndrome,...)
- Calcifying tendonitis of the shoulder
- Plantar fasciitis
- Calcaneal spurs
- Heel spur
- Iliotibial band tendonitis

For physiotherapy practices, the **Masterpuls MP100** is a very interesting choice :

"Noninvasive extracorporeal shock waves therapy (ESWT) is proving to be a successful alternative for specific indications. It is an alternative that is both medically and economically very interesting."

(Richard THIELE, MD).

In addition, shock wave therapy is also used for myofascial pain by treating trigger points. It can also be used on acupuncture points.

A special feature of the MP100 is the "vibrating massage", which will relax the muscle when applied to a contracture point.

GymnaUniphy
YOUR PHYSIO COMPANY

www.gymna-uniphy.com

info@gymna-uniphy.com